**SCIENCE
SERVICES**

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About us

Since 1976, SCIENCE SERVICES has been an excellent partner in providing systems and materials for sample preparation in microscopy and ultracentrifugation. We are trusted by scientific and industrial laboratories throughout the globe.

Since its foundation SCIENCE SERVICES works with a high degree of technical understanding when it comes to manufacturing hairpin filaments and distributing highly specialised laboratory equipment.

Today SCIENCE SERVICES offers over 20,000 products from a select group of specialist manufacturers with whom a longstanding relationship has been formed.

The core of our portfolio are products for microscopy and ultracentrifugation for material and life sciences. However, many of our products are also used in other areas.

Our team consists of a good portion of scientists as well as application specialists, technicians, engineers and great, professional staff in the office and shipping.

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▶ **Cutting Speed**

0.1 to 100mm/s cutting speed adjustable in 0.1mm/s steps throughout whole speed range.

▶ **Scan and Tilt**

A small knob on the microscope turret enabling to scan the knife edge area without turning the entire microscope even at high magnifications.

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PT PC with HD video package and view for training and documentation

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The PT PC has all the features of the PT XL and is further equipped with PC control and camera for live image acquisition and a video and snap-shot function. The PC control enables documentation and archiving of user settings and protocols. The online support using internet capabilities is included, for a direct link to specialists to assist troubleshooting.

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The high definition video package to monitor the sectioning process in real-time is displayed at the center of the PC monitor. An image capture system allows videos or still photographs to be taken at any time for data sharing in

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the lab or for collaboration as well as for publications or archiving.

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The ultramicrotome can be controlled with the touch screen monitor the digital tactile controller or with both.

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▶ **Dual Drive**

Touch screen PC and digital tactile controller.

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▶ **Video Monitoring**

The center of the touch screen displays the sectioning process in real time. Essential for user trainings and demonstrations. With recording function for acquisition of video and snap shots for documentation.

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▶ **VisuTrac**

Define and adjust cutting window during sectioning with a single mouse click.

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▶ **Digital Measuring Tool**

Define and store sample and block dimensions.

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▶ **Trimming Function**

Automated trimming into the sample block with defined depth (max. 200µm).

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▶ **Integrated Report Generator**

Direct read-out on PC to user storage via ethernet, without the need of portable data storage (e.g. USB stick).

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▶ **Online Trouble Shooting**

Skype access to RMC Service engineers.

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
PT 3D - first 1mm advance ultramicrotome

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For decades, a specimen advance of max. 200µm was the standard in ultramicrotomy. With the newly designed cutting mechanism, the PT 3D sets a new standard and is the first ultramicrotome with the ability to cut ultrathin sections up to one millimeter into a sample without resetting the advance mechanism.

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The unique 1mm specimen advance is not only beneficial for ultramicrotomy applications such as 3D reconstruction of materials, tissues, cells or cellular structures using array tomography, but also for those involving frequent trimming with the highest precision, like for the sample preparation for cryosectioning.

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Identical to the PT PC the PT 3D has a high definition video package to display and control the section process in real-time on the computer monitor. The PC software further includes all other functions, e.g. the image capture system, the digital measurement tools and more.

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► Reduced Handling

The rate of interfering with knife and sample because of resetting the advance is reduced by five times.

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► Every section is sacred

For 3D reconstruction, it is crucial to reduce the number of lost sections to a minimum. The main risk of losing a section is the re-approach to the knife edge.

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► Less Knife Damage

Most damage to a diamond knife occurs during approach of knife and sample.

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
► **RMC Specimen Holders**

Made of sturdy stainless steel with side entry for backlighting the sample in the block. Backlighting also possible for pre-trimming with the trimming block. With diameter of 10 mm, compatible with ultramicrotomes of common manufacturers.

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► **DiATOME ultra 45°, 3mm**

Standard knife for routine work. High knife stability with higher compression due to 45° angle.

DU4530	DiATOME ultra 45° 3,0mm <i>DiATOME ultra 45° 3.0mm</i>
--------	---

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► **DiATOME ultra 35°, 3mm**

Delicate knife for thinner sections and series. Smaller compression with more fragile cutting edge.

DU4530	DiATOME ultra 45° 3,0mm <i>DiATOME ultra 45° 3.0mm</i>
--------	---

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► **Vibration Isolation table by TMC**

Compact design and optimized ergonomics, compressed air connection required.

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LN Ultra - Cryo chamber for highest temperature stability

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Sectioning of both biological and industrial materials at cryogenic temperatures can be achieved with the LNUltra a liquid nitrogen-cooled cryo attachment for the PowerTome.

Capable of holding temperatures down to -180°C , the stable environment of the LNUltra is perfect for sectioning and planing a variety of samples including rubber and soft polymers, as well as biological tissues and cells, e.g., for immunolabelling using the Tokuyasu technique or CEMOVIS (Cryo Electron Microscopy of Vitreous Sections).

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Due to the low Liquid Nitrogen (LN2) consumption, a small 12-liter table-top Dewar is sufficient to provide a stable sectioning environment for around 6 hours, while saving LN2 costs.

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In combination with the PowerTome 3D, sample trimming and cutting can be achieved without resetting the specimen advance and without the need for realigning sample and knife.

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► Low LN2 Consumption

Approx. 2L/h with use of a space-saving 12L tabletop dewar.

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► High Temperature Stability

From $+20^{\circ}\text{C}$ to -180°C due to use of exclusive materials.

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► In-Chamber Specimen Freezing

Open space in the chamber for easy freezing or mounting the specimen on the pin.

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► **Easy Chamber Access**

Reduced depth of chamber facilitates sample handling and manipulation.

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► **Workstation**

Lay up tools and hands for working at the chamber without becoming cold.

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
Reduced electrostatic charging of specimens

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DCIMM3530

► **DiATOME cryo immuno 35°, 3mm for Tokuyasu**

Specially manufactured for sample preparation using the Tokuyasu technique

DiATOME cryo immuno 35° 3,0mm, für Tokuyasu

DiATOME cryo immuno 35° 3,0mm, for Tokuyasu

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DCO3530

► **DiATOME cryo dry 35°, 3mm**

for dry sectioning of material samples, e.g. polymers and films

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Diatome trim 45°, for precise trimming of sample blocks at room and cryo temperatures

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E71952-01

► **Blades**

*Rough pre-trimming of sample blocks. **Caution: Only use cleaned blades***

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Single Edged Razor Blades for Clean Rooms, Surgical Stainless Steel

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Double Edged Razor Blades, Feather®, Carbon Steel

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Dumont Style 5, straight Tweezer for handling of TEM-Grids and samples

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Dumont Style 7 curved Tweezer for handling of TEM-Grids and samples

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Array Tomography in Electron Microscopy

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
Array Tomography in Electron Microscopy is one way to acquire data for the 3D reconstruction of structures. Thereby, serial sections are produced and collected on a substrate for subsequent imaging. Depending on the size or volume of a targeting region more or less sections are required for reconstruction.

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Array Tomography allows imaging of the same section at different pixel resolutions. Images taken of low resolution (700 nm/pixel), intermediate resolution (30 nm/pixel), and high resolution (3 nm/pixel).

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In comparison to other 3D-techniques, such as Serial Block-Face Scanning EM (SBFSEM) and Focused Ion Beam (FIB), Array Tomography is non-destructive. Collected sections can be stored, re-imaged at difference regions and resolutions, re-stained with heavy metals to enhance contrast or with fluorescent markers for correlative light and electron microscopy (CLEM).

Array Tomography is used in various research fields not only for reconstruction of biological structures but also for complex materials.

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► **Non-Destructive**

Sections are stored on substrate and can be kept for years.

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► **Sample Library**

Stored sections can be re-evaluated at different locations and resolutions.

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


► **Save Scan Time and Data Storage**



No need for highest resolution or whole area scans due to re-imaging at various resolutions and regions.

► Large Volume Reconstruction with ATUMtome



▪ (a) Image of cortical column neurons in the cerebral cortex. (b,c) About $1000\mu\text{m}^3$ around an apical dendrite of one pyramidal cell. This region contains around 675 synapses, 530 axons, and 90 dendrites with highly non-random connectivity.

▪ Images of tissue sections reconstructed through 3D imaging and array tomography.


*Images supplied courtesy:
Jeff Lichtman, Harvard University,
Boston, USA*

► Small Volume Reconstruction with ASH

▪ Root tip of *Arabidopsis thaliana*. (a) LM and EM image of section ribbons on ITO substrate. (b) Reconstruction of root tip vacuoles. Large vacuole in yellow wire mesh, small vacuoles in light blue. Volume rendering performed with UCSF Chimera.

*Images supplied courtesy:
Irene Wacker and Rasmus Schröder, BioQuant, Heidelberg, Germany.*

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ATUMtome – A fully automated system for large volume reconstructions

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In classical ultramicrotomy, serial sections are cut and float on the water surface in the diamond knife trough and are picked up manually.

The ATUMtome is a unique system for automated collection of sections on a continuous tape. Thereby, the tape moves through the knife trough and automatically collects the sections from the water surface in sequential order. After section collection the tape is cut in strips and is mounted on substrates such as silicon wafers or glass slides. Data acquisition is done with a scanning electron microscope.

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The ATUMtome in conjunction with the Powertome 3D can produce many thousands of sections uninterrupted thanks to the 1mm advance of the PT3D. This means no loss of sections which is crucial for array tomography.

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► **Fully Automated**

The tape collection and the cutting process are synchronized. Tape and cutting speed can be adjusted individually.

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► **Unattended Section Collection**

Fully automated section collection process. Once set and stable, section pick-up by tape without the need of human intervention.

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► **Easy Handling**

Intuitive and simple section collection mechanism. No long training period is required.

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► **Keeps Microscope Free**

Sectioning process does not interfere with image acquisition. In SBFSEM or FIB sectioning is executed in the microscope.

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Controlled Environment Considerations

Stabile Temperatur: ± 1°C /h <i>Stable temperature: ± 1°C /h</i>	Abschirmung gegen Windströmungen, Lüftungsanlagen <i>Shield from wind currents, ventilation systems</i>
Rel. Luftfeuchtigkeit: ~50% <i>Rel. Humidity: ~50%</i>	Verhindern der Probenblock-Aufladung beim Schneiden (DiATOME Ionisator) <i>Prevent sample block charging during sectioning (DiATOME Ionizer)</i>

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DiATOME trim 90, for orthogonal trimming for serial sectioning

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DiATOME ultra 35° 4.0mm, for sectioning of large sample blocks

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ASH2 – The Advanced Substrate Holder of the 2nd Generation

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The ASH2 is a multi-axis manipulation tool designed to precisely position substrates in a knife trough and is capable of collecting tens to hundreds of serial sections for array tomography and correlative 3D workflows.

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A substrate is clamped to the ASH2 and placed in the water of the knife trough. After sectioning, the water is slowly removed and the section ribbons adhere to the substrate.

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The easy setup, alignment and handling of the ASH2 allows an efficient collection of sections on a desired substrate even for small section series. Subsequent data acquisition is executed at the light (LM) or scanning electron microscope (SEM).

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The main application of the ASH2 is the take-up of a few hundred sections for reconstruction of small sample volumes.

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► **Fast Set Up**

Quickly mounts onto the knife stage, converting a standard ultramicrotome to an array tomography system.

► **No Mounting Rail**

No need for a mounting rail and special table like for the precursor ASH-100.

► **Fits Any Ultramicrotome**

Integrates with current RMC and Leica ultramicrotomes including UC6 and UC7.

► **No Limitation to Knives**

Compatible with commercially available diamond knives for room temperature. Usage with a variety of standard diamond knives as well as knives with larger troughs for large substrates.

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A simple solution for manual collecting of section ribbons is the ultra ats knife. This DiATOME diamond knife is characterized by a particularly large water trough that can be mounted with a substrate. The ribbons are manually maneuvered over the substrate. By lowering the water level through a hole in the bottom, the sections adhere to the substrate and are ready for further processing.

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RMC Boeckeler MT990 - Rotary Microtome with CR1000 cryo-sectioning system

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The RMC MT990 series of microtomes are made for materials science as well as histology, for room temperature and cryo applications. The MT990 rotary motorized microtome provides a reproducible, even cutting speed for the hardest of samples. The MT990 accepts all common knife types in microtomy, including glass, diamond as well as disposable blades, solid steel or tungsten carbide knives.

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In conjunction with the CR1000, the MT990 becomes the world's only cryo microtome. It is the ideal tool for cryo-sectioning samples for optical and FTIR transmission microscopy as well as cryo-planing sample surfaces for AFM and SEM. Large area rubber, polymer and wet biological samples can be cut at temperatures down to -160°C with diamond and glass knives.

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The CR1000 includes most of the features and technologies used in cryo-ultramicrotomy:

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- Accepts all common knife types
- Adjustable LED overhead- and backlight
- Knife rotation $\pm 20^\circ$
- Liquid Nitrogen Dewar with 12l capacity for up to 6 hours operation
- Independent temperature controls of knife and specimen
- Temperature range +40 down to -160°C
- High temperature stability of $\pm 0.1^\circ\text{C}$

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► Cutting Stroke

The MT990 has a cutting stroke of 64mm. With CR1000 cryo-attachment the cutting stroke is reduced to 40mm. The cutting stroke has to be reduced at the factory.

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Knives and Blades

The sectioning process is not complete without knives and we provide cutting edges for all types of EM and LM applications. Tungsten carbide knives are popular for a number of hard sample applications in materials science and diamond knives are the preferred tool for delicate applications. For many, the glass knife is preferred for teaching, trimming and for those unknown samples which may have inclusions that can instantly damage your diamond edge. Disposable blades are standard for paraffin embedded and other soft specimens.

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
Tungsten carbide blade, box of 3, 38mm wide x 12.7mm high x 1mm thick. For semi-thin sectioning with the CR 1000.

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Tungsten carbide knife, Profile D, 16cm, with case. For use on MT990 Type N knife holder

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MS200 Disposable universal blade, 50 pcs, MICROS Stainless steel blade with plasma coating for daily routine, suitable for all types of tissue, especially for ribbon sectioning and thin cuts. Angle: 35°

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MS100 Disposable blade, 50 pcs MICROS Stainless steel blade with teflon coating for daily routine, especially for hard and difficult tissue, suitable for cryostat cutting. Angle: 35°

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
MS33 Disposable blade, 50 pcs, MICROS Stainless steel blade for very thin sectioning of all types of tissue, also suitable for skin biopsy. Angle: 30°

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The vibrating microtomes of the 7000- and 5100-series have been used in the fields of electrophysiology, neuropathology and immunohistochemistry for more than 30 years. The vibrating blade allows fixed tissue sectioning and the preparation of acute tissue slices. The use of fresh tissues results in less artifacts or changes of morphology. The oscillation reduces compression, cell destruction or other incisional damage.

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7000smz-2 – high precision sectioning for optimal tissue slicing

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The advanced high precision vibrating microtome 7000smz-2 provides optimal tissue sectioning for visual patch-clamp recording or high resolution imaging. With a minimal z-axis deflection of less than 1µm the vibratome delivers healthy viable sections of all biological tissues. The 7000smz-2 is supplied with its own Z-axis calibration verifier.

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In addition to user-friendly operation, the versatile user interface offers a range of adjustable slicing parameters such as cutting head advance speed, blade oscillation frequency or section thickness.

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- ▶ **Precise Section Thickness increments**
Adjustable section thickness step size of 1µm.
- ▶ **Customizable Blade Holder Angle**
Blade holder angle individually adjustable according to user requirements.

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▶ **Included Starter Kit**

Starter kit with all required consumables (blades,..) .

▶ **Operation Modes**

Manual, semi or fully automatic operation.

▶ **Parameters & User Profiles**

Saving and adjusting of slicing parameters with 8 customizable user profiles.

▶ **Minimal Z-axis Deflection**

Less than 1 µm at all speeds and amplitudes with included z-axis calibration unit.

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▶ **Tissue Bath**

Temperature controlled bath with controllable temperature range of approx. +50°C to +10°C (± 1°C accuracy)

E7611-A	Temperaturgesteuertes Standard-Gewebebad für 7000-Modelle <i>Temperature Controlled Standard Tissue Bath for 7000-series</i>
---------	---

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
The 5100mz and 5100mz-plus are designed for routine techniques such as histology, organotypic slice cultures and low resolution imaging. The 5100mz-plus includes a z-axis calibration unit, which reduces damage to the slice surface and provides longer viable slices required for e.g. electrophysiological field recordings or ISME's.

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The 5100mz-series share many features with the high precision vibrating microtome 7000smz-2 such as the versatile and operator-friendly user interface.

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► **Minimal Z-axis Error of Less than $8\mu\text{m}$**

The optional calibration unit reduces the z-axis blade deflection to $2\mu\text{m}$ (5100mz-plus incl.).

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► **Manual or Semi-automatic Operation**

Semi-automated mode automatically remembers start and finish position of the slice.

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► **Sterile Use in Culture Hood**

Remote controls allowing for operation outside the hood.

Bestellinformationen / Order Info

Artikelnummer / Product Code	Bezeichnung / Description
E5100mz	5100mz Vibratom 5100mz Vibrating Microtome
E5100mz-Plus	5100mz-plus Vibratom 5100mz-plus Vibrating Microtome

Zubehör 5100mz Serie
Accessories 5100mz series

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▶ **LED Light & Magnifying Glass**
Integratable Cold Light Source & Magnifier - approx. 1.5x magnification

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▶ **Upgrade Kit for 5100mz**
Optional kit to upgrade the 5100mz to the 5100mz-Plus; optimally calibrate the Z-axis deflection of the blade to 2µm

E5100-1-45

Upgrade Kit für 5100mz
Upgrade Kit for 5100mz

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
Vergleichstabelle - Vibratome

Comparative Vibrating Microtomes

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The GKM2 glass knife maker produces the most reproducible straight edged glass knives available. With its balanced break methodology coupled with a pressure detector readout for precise reproducibility of breaking pressure, the GKM2 is the favored knifemaker for many laboratories.

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► Produce Optimum Knives

Produce knife edges 6.4mm or 8mm wide with different angle from 35° to 55°.

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► Easy Positioning

Accurately and quickly position glass strips with flip down positioning stops. Unique positioning micrometers to accurately position the glass square for final score.

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► Balanced Break

The GKM2 is specifically designed to incorporate the „balanced break“ method to produce optimum knife edges. Balanced breaks require equal lengths of glass on each side of the score, along with a mechanism that applies equal load to each side.

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► Controlled Breaking Force

To achieve consistent, controlled breaks, two aspects of the breaking process must be measurably controlled: position and load. The load applied to the glass is controlled through the measurement and constant feedback displayed on the digital readout as the clamping handle and breaking arm are adjusted.

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With these unique features it is possible to apply reproducible force from knife to knife with the GKM2 – to break the glass as slow as you like.

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
Note: Cleaning The Glass Strips

The glass strips should be carefully cleaned with a good quality mild detergent such as used for cleaning fine crystal glassware. Wear rubber gloves and be careful handling the glass to avoid cutting your hands with the sharp edges. Use a soft sponge to clean off any residue. Rinse with distilled water and dry with a good quality lint-free laboratory tissue.

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Messer	Messertyp	Messerwinkel	Grösse (mm)	Schnittdickenbereich (nm)	neu #	Umtausch #	Nachschliff #
ultra							
	ultra 45°	45°	1,5	30 – 200	DU4515	TDU4515	RDU4515
			2,0		DU4520	TDU4520	RDU4520
			2,5		DU4525	TDU4525	RDU4525
			3,0		DU4530	TDU4530	RDU4530
			3,5		DU4535	TDU4535	RDU4535
			4,0		DU4540	TDU4540	RDU4540
	ultra 35°	35°	1,5	30 – 200	DU3515	TDU3515	RDU3515
			2,0		DU3520	TDU3520	RDU3520
			2,5		DU3525	TDU3525	RDU3525
			3,0		DU3530	TDU3530	RDU3530
			3,5		DU3535	TDU3535	RDU3535
			4,0		DU3540	TDU3540	RDU3540
	ultra jumbo 45°	45°	2,0	50 – 500	DUJ4520	TDUJ4520	RDUJ4520
	ultra jumbo 35°	35°	3,0		DUJ4530	TDUJ4530	RDUJ4530
			2,0		DUJ3520	TDUJ3520	RDUJ3520
			3,0		DUJ3530	TDUJ3530	RDUJ3530
	ultra semi	35°	3,0	50 – 500	DU3530-semi	TDU3530-semi	RDU3530-semi
	ultra AFM	35°	2,0	15 – 100	DU3520-AFM	TDU3520-AFM	RDU3520-AFM
			3,0		DU3530-AFM	TDU3530-AFM	RDU3530-AFM
	ultra sonic (mit Steuereinheit)*	35°	3,0	15 – 100	DUS3530	-	-
	ultra sonic (nur Messer)*				DUS3530	-	RBUS3530
cryo							
	cryo 35° (trocken)	35°	1,5	30-300	DCO3515	TDCO3515	RDCO3515
			2,0		DCO3520	TDCO3520	RDCO3520
			2,5		DCO3525	TDCO3525	RDCO3525
			3,0		DCO3530	TDCO3530	RDCO3530
			3,5		DCO3535	TDCO3535	RDCO3535
			4,0		DCO3540	TDCO3540	RDCO3540
	cryo 45° (trocken)	45°	1,5	30-300	DCO4515	TDCO4515	RDCO4515
			2,0		DCO4520	TDCO4520	RDCO4520
			2,5		DCO4525	TDCO4525	RDCO4525
			3,0		DCO4530	TDCO4530	RDCO4530
			3,5		DCO4535	TDCO4535	RDCO4535
			4,0		DCO4540	TDCO4540	RDCO4540
	cryo 35° (nass)	35°	1,5	30-300	DCM3515	TDCM3515	RDCM3515
			2,0		DCM3520	TDCM3520	RDCM3520
			2,5		DCM3525	TDCM3525	RDCM3525
			3,0		DCM3530	TDCM3530	RDCM3530
			3,5		DCM3535	TDCM3535	RDCM3535
			4,0		DCM3540	TDCM3540	RDCM3540

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The PT FL - A PowerTome with fluorescence location system

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The new RMC Boeckeler PTFL - developed from a collaboration between the Francis Crick Institute and RMC - is an ultramicrotome-mounted microscope, capable of imaging fluorescent features in a resin block during trimming and sectioning.

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With the PTFL, resin embedded samples can be quickly evaluated and the ROI (region of interest) targeted. The operator can track the ROI during trimming until the fluorescent features reach the block surface.

The researcher can then collect ultrathin sections of only the ROI for TEM (on grids) and SEM (e.g. silicon).

For Serial Block Face SEM and Focused Ion Beam SEM, the trimmed blockface containing the ROI can be removed and mounted for volume EM.

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► Save Time

On-the-fly targeting of fluorescent ROIs.

► Save Costs

Reducing user time, microscope time and data storage.

► Save Storage Space

Data is only acquired from the ROIs, with big savings in data storage and handling requirements.

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Tech@ScienceServices.de

Any further Questions?

If you have any further questions regarding a product or the field of microtomy, please do not hesitate to contact us. We also offer on-site demonstrations.

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Science Services GmbH