

Endless Possibilities ...

Köhler Alignment for the Optical Light Microscope (OLM)

The brain is an excellent image processor and can ignore image defects in real time but when viewing a captured image, that image stands alone and the brain then sees any defect. Performing a Köhler alignment will ensure that most common defects, such as uneven illumination, will not be present in the final image. However, the process is all too often an ignored concept in OLM and should be performed before every viewing session, especially when capturing images. NOTE: The most basic OLMs may not have the optical elements that are needed in order to perform the alignment.

A properly executed Köhler alignment produces an extremely even illumination of the sample and ensures that an image of the illumination source is not visible in the resulting images.



Protocols:

Alignment

Köhler

Directions

1. Select low power objective (10x) and place slide on stage.



2. Set IPD (inter-pupillary distance) © 2018 E



3. Focus image by: 3a. Bring stage up close to objective lens. 3b. Looking through oculars move specimen stage DOWN until a focused image is obtained. **Out of Focus**

irsch noteS



In Focus





- opy Sciences 4. Adjust dioptric eyepiece for proper focus:

 - scopy Sciences
 - © 2018 1



5. Open iris or condenser aperture wide open.



7. Adjust condenser height until the edges of the field diaphragm come into sharp focus. .∪ s © 2018 Electron Micros 2018 Electr

© 2018 Electron Microscopy Sciences

© 2018 Electron Microscopy Sciences

Close field 6. diaphragm all the way.



(Micros

© 2018 Electron **Image with Field Diaphragm Closed**





8. Center illumination using the condenser centering screws. Open the field diaphragm some to be sure of position.









 Close the iris or condenser diaphragm down to the point where there is a marked increase in the image sharpness or contrast (Image B).
NOTE: This should match the Numerical Aperture on the lens.

1111

9. Open field diaphragm until it is just outside the field of view. (Image A).





3

8+ 0

11. **NOTE:** Steps 6-9 should be done for each objective used, but always start at the lower magnifications first.

R (resolution) = $0.612 \times \lambda$ NA

NA (numerical aperture) = n (sin α) n= index of refraction

Electron Microscopy Sciences | 1560 Industry Road | Hatfield, PA 19440 P 215-412-8400 | F 215-412-8450 | info@emsdiasum.com