



Endless Possibilities ...

Biological Processing Steps for TEM

Chemical processing of biological tissues for microscopic examination has evolved to keep pace with increased levels of detail seen with newer technologies. From simple formalin fixation, paraffin embedment and 0.5 μ m sections for OLM into a multichemical, epoxy resin embedment and 60 nm sections for TEM. Even this well-established TEM process has evolved further with the use of microwaves to significantly decrease the time and potential processing artifacts. This generalized procedure gives the acceptable changes, concentrations and times for both bench top and microwave processing.

Setup requirements

- Program microwave for desired processing times.
- Get ORGANIZED!! Have equipment and solutions ready.
- Samples should be placed on a rotator when using Benchtop processing.

All of the following steps can be carried out in the 1.7 ml microfuge tubes, scintillation vials, or Petri dishes. **NOTE:** Acetone or propylene oxide (PO) cannot be used for dehydration if plastic Petri dishes are used.

STEP	- pol Se	S	TEMP	Microwave	Lynx II/Benchtop
Nicro	Initial fixation (Karnov	vsky's)	37º C	2:30 min	2 hr. cient
2.	Buffer rinse 3 changes		37º C	60 sec. ea.	10 min. ea.
3.	2-4% OsO_4 in DI wat (Sometimes 2% Pota	er ssium Permang	37º C anate in DI is us	2:30 min ed for plants and bacteria.)	2 hr.
4.	Water rinse 3 chang	ges	37° C	60 sec. ea.	10 min. ea.
NOTE	: If using LR White	e, Acetone or	PO can not l	be used, only ETOH!	ciences
5. 29	Dehydration (Using either:	50% 70%	45° C 45° C ∽	60 sec. 60 sec.	10 min. 10 min.
	ETOH, Acetone or Acetonitrile)	80% 90%	45° C 45° C	60 sec. 60 sec.	10 min. 10 min.
- citi	2 changes	100%	45° C	Elect 60 sec. ea	Micros ⁰ 10 min. ea
NOTE	: Separate the SE	M for CPD or	HMDS from t	he TEM samples at this	stage if necessary.
(ETOH	: Acetone OR	1:1 100%	45° C 45° C	60 sec.	5 min.
ЕТОН	: PO)	internet		© ² goo.	Science
6.	Infiltration - ETOH, Acetone, or PC		Resin	cC	0,67
	Plan	t 3:1	50° C	15 min.	30 min.
	10 ^{SU}	2:1	50° C	5 min.	1 hr.
	Mici	1:1	\50° C	15 min.	1 hr.
	100% resin	2 changes	⊂ 50° C	15 min. ea.	1 hr. ea
78E	Embed in capsules and polymer Or store in freezer in 100% resir		over night at 70° il time for Embe	dment. EMS Catalog S	Supplies listed on back
	-018 E.				i oction.

